



## Segregation analysis offers a mechanism for variant reclassification in a small subset of cases but is especially powerful in classifying deleterious mutations

Julie Eggington, Jennifer Saam, Karla Bowles, Kelsey Moyes, Susan Manley, Lisa Esterling, Scott Sizemore, Christopher Arnell, Brandon Robinson, Eric Rosenthal, Jeremiah Bennett, Mark McKellar, Benjamin Roa and Richard Wenstrup\*

Poster presented at the European Society of Human Genetics (ESHG) Congress 2013


### Introduction

Accurate classification of variants in regards to their clinical significance is a critical challenge associated with gene sequencing tests, with this challenge expected to increase dramatically as next generation sequencing technologies are more widely used. Thus, it is important to evaluate the effectiveness of various strategies for variant reclassification.


At Myriad Genetic Laboratories Inc., we utilize multiple lines of evidence to evaluate and reclassify variants of uncertain significance (VUS), and new classification methods are continually being assessed by Myriad scientists. Here we describe the use of segregation analysis for the reclassification of variants in the genes associated with Hereditary Breast and Ovarian Cancer syndrome (HBOC). We report on the participation rate of Myriad's family testing for segregation analysis through its Variant Classification Program.

We also show that for the discovery of benign variants, segregation analysis is a comparatively weak method compared to other methodologies developed at Myriad. However, for the identification of truly deleterious mutations, we have demonstrated that segregation analysis is a powerful method and compliments the other variant reclassification methods employed at Myriad.

## Reclassification Methods in Myriad's Variant Classification Program



Frequency of use by Reclassification Method  
Oct 2011 – Nov 2012




Percentage of Reclassification Upgrades and Downgrades Achieved by Each of Myriad's Reclassification Methods  
Oct 2011 – Nov 2012




Relative Effectiveness of Reclassification Methods by Proband Count  
Oct 2011 – Nov 2012

Reclassification Method	Average Number of Proband per Downgrade & (minimum)**	Average Number of Proband per Upgrade & (minimum)**
<i>in trans</i>	17.5 (1 min)	NA
Mutation Co-Occurrence	31.2 (2 min)	NA
Segregation	41.8 (9 min)	32 (13 min)
History Weighting Algorithm	19.4 (6 min)	33.7 (28 min)
Evolutionary Conservation	15.4 (NA)	NA
Literature Review	7 (NA)	12.6 (NA)
Population Frequency	8.5 (NA)	NA

\*\* These represent the real life minimum data requirements necessary for reclassification in the Oct 2011 – Nov 2012 time period. They do not necessarily represent theoretical minimums.



Decline in Rate of BRCA1/2 Variants of Uncertain Significance found in USA Patients



# Segregation analysis offers a mechanism for variant reclassification in a small subset of cases but is especially powerful in classifying deleterious mutations

## Conclusions:

The majority of *BRCA1/2* VUSs are discovered to be benign through a variety of methods, with History Weighting Algorithm being the most robust method. Segregation analysis shows particular power in identifying deleterious variants rather than benign variants. Considering laboratories and community research centers have finite resources, these data therefore suggest that the tailoring of family analysis to specific families with higher likelihoods of having a deleterious mutation may be the most productive use of resources.

## References

1. Hall MJ et al. *Cancer*. 2009. 115(10):2222-2233.
2. Frank TS et al. *J Clin Oncol*. 2002. 20(6):1480-1490.
3. Tavtigian SV et al. *Familial Cancer*. 2006. 5(1):77-88.
4. Easton DF et al. *Am J Hum Genet*. 2007. 81(5):873-883.

Concurrent ESHG 2013 Poster number P11.047: "A Clinical History Weighting Algorithm Accurately Classifies *BRCA1* and *BRCA2* Variants, Bowles et al".

## BRACAnalysis®

A test for Hereditary Breast and Ovarian Cancer (HBOC) syndrome

A predictive medicine product for hereditary breast and ovarian cancer.

BRACAnalysis® testing assesses a woman's risk of developing hereditary breast or ovarian cancer based on detection of mutations in the *BRCA1* and *BRCA2* genes. This test has become the standard of care in identification of individuals with hereditary breast and ovarian cancer.

\*All authors are employees of Myriad Genetics, Inc., and Myriad Genetic Laboratories, Inc. and receive salary and stock options as compensation.



MYRIAD®

Myriad Genetics GmbH  
Leutschenbachstrasse 95  
8050 Zurich  
Switzerland

[www.myriad.com](http://www.myriad.com)